
technical manual

MANUALE TECNICO
TECHNICAL MANUAL
MANUEL TECHNIQUE

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**ABBREVIATIONS,
TEMPERATURE AND CHEMICAL
RESISTANCE OF PLASTICS**

The table below lists commonly used abbreviations for plastics. This list covers plastics commonly employed in the manufacture of plastic laboratory ware.

Temperatures appearing in parentheses: limits tolerated for intervals only.

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"CLEAN" AND ECOLOGICAL
PLASTICS

Thanks to technological innovation and progressive adaptation to the most elevated eco-compatibility standards, Kartell plastic raw material are of excellent quality. In fact the used plastic materials grant wide compatibilities; hereunder you can find some standards conformities related to their characteristics.

- Foodstuff contact suitable (National and European Directive)
- Foodstuff contact suitable (FDA Directive)
- PTFE: foodstuff contact suitable (FDA CFR TITLE 177.1550)
- Absence or limitation according the International Directives of: heavy metals, biphenyls and phtalates
- RoHS Directives Conformity
- BSE and/or TSE Directives Conformity
- Atex Directives Conformity

Kartell Labware promotes the research in the field of the plastic materials, through the direct contact with the most important worldwide raw materials manufacturers and the research of innovative materials. We would like to underline that Kartell Labware materials, if not chemically contaminated, are totally recyclable

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STERILISING PLASTIC LABORATORY WARE

Before sterilising any items of plastic laboratory ware, verify that no contamination or residues are present. Their presence could destroy plastics during sterilisation or autoclaving.

Observe the tolerated temperature ranges of plastic when autoclaving plastic laboratory ware. Remove any stoppers, fittings, or caps from plastic laboratory ware prior to autoclaving. Plastic vessels should be autoclaved separately from their closures and other fitting. Autoclaving vessels with their closures in place can lead to deformation and destruction of the vessels. All statements are advisory only, and imply no liability on the part of Kartell.

All statements relating to the resistances of plastic laboratory ware to high temperatures, chemicals, and to sterilisation and cleaning procedures have been cautiously formulated, based on statements of raw material manufacturers, on statements appearing in the literature, plus experience gained in actual practice.

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CLEANING PLASTIC LABORATORY WARE

All polyolefins, such as LDPE, HDPE, PP and PMP (TPX®), as well as the fluorinated hydrocarbons PTFE, PFA, FEP, ETFE and E-CTFE, have wettable surfaces that are both highly resistant to high temperatures and chemical attack and easy to clean. Slight contamination can be removed using a chemically neutral (pH 7) cleaning agent. Heavy contamination can be removed using an alkaline (pH up to 12) cleaning agent.

Use only chemically neutral (pH7) cleaning agents on polycarbonate (PC) or Polystyrene (PS).

Laboratory ware used in trace analyses should be cleaned in a 1-N hydrochloric acid (HCl) solution for periods of not more than 6 hours and then rinsed in distilled water in order to preclude contamination by cations or anions.

Never use scouring powders or abrasive sponges when cleaning plastic laboratory ware. Use no alkaline cleaning agents on polycarbonate (PC) laboratory ware.

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